

## NuboMed S1000 Height Adjustable Cart



NuboMed

- **Material:** Constructed from a combination of high quality PC & ABS and aluminium alloy
- **Integrated Design:** Handle and work surface seamlessly integrated from a single piece of anti-microbial PC & ABS material, ensuring a gap-free structure inclusive of a spill-proof groove design
- **Ergonomic and Effortless:** A simple, non-powered design for seamless integration during patient rounds
- **Height Adjustable:** Pneumatic adjustment with a range of 0-300mm
- **Base:** Designed to prevent tripping and made from a single-piece of aluminium alloy with an antibacterial powder coating finish
- **Optional Accessories:** Enhanced functionality through a range of compatible accessories and modifications
- **Castors:** Equipped with 2x 4' lockable and 2x 4' tracking castors
- **Ideal For:** Doctor's Round/Training/Conferencing, Daily Nursing



### Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height
WAT908005	NuboMed S1000 Height Adjustable Cart	576 mm	765 mm	1563 mm