


## Lateral Chest Support Assembly

- A great range of accessories to fit the Aspire Manual Wheelchair range to improve functionality and user experience
- Click [here](#) to ensure this accessory fits your specific wheelchair model



### Warranty Information

- 1 Year

 Contract approved

**Compliance**

AS3696-1:2008



Code	Display Name	RRP	Box Quantity
WCA848430	Lateral Chest Support Assembly - Left	127	1
WCA848450	Rehab RS/RX Wheelchair Lateral Chest Support Assembly - Left	203	1
WCA848440	Lateral Chest Support Assembly - Right	127	1
WCA848460	Rehab RS/RX Wheelchair Lateral Chest Support Assembly - Right	203	1