

Meiko Bottle Rack



- Drying and storage rack for urine bottles
- Stainless steel
- Rack at the bottom
- Drip tray is sold separately, size is compatible with storage rack

 Contract approved

Support @ Home

Function: Collecting urine and faeces

Category: Prescribed

Click [here](#) for more information



Code	Display Name	Total Length	Total Width	Total Height	Capacity
UPA400210	Meiko Bottle Rack - TU70 - 11 Urine Bottles	350 mm	700 mm	650 mm	11 Urine Bottles