

ROHO - Hybrid Elite Sensor Ready



permobil

ROHO

The HYBRID ELITE SR Single Compartment cushion combines air and foam in one cushion and is designed for individuals whose primary need is stability. When sitting on the cushion, your pelvis is held and protected with the air insert, and the thighs are aligned with the contoured foam leg channels. Using the available Smart Check device, you can inflate your cushion to the proper level and check it whenever you like with the simple push of a button.

- **The stability of foam. The protection of air. All in one**
- Stabilises pelvis and thighs for ideal positioning
- Easy to set up and can be used with the optional Smart Check device for set up and monitoring
- Adjustable to shape and weight in the pelvic loading area
- Adjustable to changes in individual's position in the pelvic loading area over time
- Hydrostatic distribution of load across the pelvic loading area
- Provides up to approximately 70 mm of immersion
- Air-filled cushion insert can be easily cleaned and disinfected
- Cover is machine washable and can be disinfected, foam base may be wiped with damp cloth
- Neoprene cushion insert can be patched or repaired to extend use
- Contoured foam base, foam insert, and air-filled insert are fluid resistant
- Foam insert is contoured to provide support to the legs
- Smart Check purchased separately (MSP500488)

Warranty Information

- 2 Years

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_181006190_0105_1_2

Pressure Care Cushion –

Composite:05_429_0103_1_2

Click [here](#) for more information

Compliance

EN-ISO 10993-10 (1999-11) / EN 1021-1 and EN 1021-2

Accessories

Click on the cart icon to display the range of optional accessories

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Pressure Care Rating
PTC594860SR	ROHO Hybrid Elite SR - Single Compartment - 14 x 16in	14 in	16 in	4.25 in	225 kg	High Risk
PTC594899SR	ROHO Hybrid Elite SR - Single Compartment - 15 x 15in	15 in	15 in	4.25 in	225 kg	High Risk
PTC594900SR	ROHO Hybrid Elite SR - Single Compartment - 15 x 16in	15 in	16 in	4.25 in	225 kg	High Risk
PTC594901SR	ROHO Hybrid Elite SR - Single Compartment - 15 x 17in	15 in	17 in	4.25 in	225 kg	High Risk
PTC594909SR	ROHO Hybrid Elite SR - Single Compartment - 16 x 15in	16 in	15 in	4.25 in	225 kg	High Risk
PTC594910SR	ROHO Hybrid Elite SR - Single Compartment - 16 x 16in	16 in	16 in	4.25 in	225 kg	High Risk
PTC594911SR	ROHO Hybrid Elite SR - Single Compartment - 16 x 17in	16 in	17 in	4.25 in	225 kg	High Risk
PTC594920SR	ROHO Hybrid Elite SR - Single Compartment - 16 x 18in	16 in	18 in	4.25 in	225 kg	High Risk
PTC594921SR	ROHO Hybrid Elite SR - Single Compartment - 16 x 20in	16 in	20 in	4.25 in	225 kg	High Risk
PTC594929SR	ROHO Hybrid Elite SR - Single Compartment - 17 x 16in	17 in	16 in	4.25 in	225 kg	High Risk
PTC594930SR	ROHO Hybrid Elite SR - Single Compartment - 17 x 17in	17 in	17 in	4.25 in	225 kg	High Risk
PTC594931SR	ROHO Hybrid Elite SR - Single Compartment - 17 x 18in	17 in	18 in	4.25 in	225 kg	High Risk
PTC594935SR	ROHO Hybrid Elite SR - Single Compartment - 18 x 16in	18 in	16 in	4.25 in	225 kg	High Risk
PTC594939SR	ROHO Hybrid Elite SR - Single Compartment - 18 x 17in	18 in	17 in	4.25 in	225 kg	High Risk
PTC594940SR	ROHO Hybrid Elite SR - Single Compartment - 18 x 18in	18 in	18 in	4.25 in	225 kg	High Risk
PTC594945SR	ROHO Hybrid Elite SR - Single Compartment - 18 x 20in	18 in	20 in	4.25 in	225 kg	High Risk
PTC594950SR	ROHO Hybrid Elite SR - Single Compartment - 20 x 18in	20 in	18 in	4.25 in	225 kg	High Risk

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Pressure Care Rating
PTC594951SR	ROHO Hybrid Elite SR - Single Compartment - 20 x 20in	20 in	20 in	4.25 in	225 kg	High Risk
PTC594946SR	ROHO Hybrid Elite SR - Single Compartment - 22 x 18in	22 in	18 in	4.25 in	225 kg	High Risk
PTC594947SR	ROHO Hybrid Elite SR - Single Compartment - 22 x 20in	22 in	20 in	4.25 in	225 kg	High Risk
PTC594948SR	ROHO Hybrid Elite SR - Single Compartment - 24 x 18in	24 in	18 in	4.25 in	225 kg	High Risk
PTC594949SR	ROHO Hybrid Elite SR - Single Compartment - 24 x 20in	24 in	20 in	4.25 in	225 kg	High Risk
PTC594845SR	ROHO Hybrid Elite SR - Single Compartment - 14 x 14in	14 in	14 in	4.25 in	225 kg	High Risk
PTC594850SR	ROHO Hybrid Elite SR - Single Compartment - 14 x 15in	14 in	15 in	4.25 in	225 kg	High Risk