

Paediatric Cot



- Paediatric Cot features a stainless steel or epoxy coated frame
- Mechanisms to operate drop sides located outside of the cot at mattress base for added safety
- Standard foam 50 mm mattress
- Four castors - two lockable for stability
- Can be tilted at Trendelenburg angle, up to 12 degrees above the horizontal
- Meets AS 2130



Compliance

No applicable compliance standard

Warranty Information

- 3 Years

Code	Display Name	Total Length	Total Width
WAT830500	Paediatric Cot - White	1440 mm	760 mm
WAT830100	Paediatric Cot - Double Latch - Stainless Steel	-	-
WAT830700	Paediatric Cot Economy Mattress	1380 mm	700 mm
WAT830600	Paediatric Cot IV Pole	-	-
WAT830900	Paediatric Cot Oxygen Bottle Carrier	-	-
WAT831000	Paediatric Cot Safety Enclosure	-	-
WAT831100	Paediatric Cot Safety Enclosure - Epoxy Coated	-	-

Code	Display Name	Total Length	Total Width
WAT764200	Cot Paediatric Mobile S/S	-	-