

Digital TENS Unit with Electrodes & Leads



A dual channel digital TENS with memory function that can store the treatment time used by the patients. The therapist can then check the units to make sure the unit is used as suggested. This device has a large liquid crystal display and press buttons, which is a great design for people who have difficulty reading the small inscription on the panel or turning the small dials.

- Channel: Dual, isolated between channels
- Pulse Amplitude: 0-80mA peak to peak into 500 ohm load
- Pulse Rate: Adjustable, from 1Hz ~ 150Hz, 1Hz/step
- Pulse Width: Adjustable, from 30 uS ~ 260 uS, 10 uS/step
- Mode: B, N, M, S1, S2
- Timer: Selectable, 15, 30, 60 Minutes
- Wave Form: Asymmetrical Bi-Phasic Square Pulse
- Voltage: 0-40 V peak into 500 ohm load
- 9 preset programs
- Supplied with the following accessories:
 - 1 x Carry Case
 - 2 x TENS Leads
 - 1 x Pack of AllCare Self-Adhesive Electrodes (4 Electrodes/pack)
 - 1 x Operating Manual
 - 1 x 9V Battery

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Not on S@H

Category: Not on S@H

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Total Height	Unit Weight
FPP315210	TENS Unit - Digital w/Electrodes & Leads	101 mm	68 mm	215 mm	155 g

