

OPPO 2321 ELASTIC KNEE SUPPORT



- 3D energising pattern stimulates sensory inputs to facilitate joint stabilisation during movement
- Elasticated zones offer freedom of movement with ideal coverage
- Comfortable edge design holds firm and provides a snug fit without being too tight
- Seamless design enhances all-around comfort and compliance Application:
 - Muscle Strain
 - Knee Sprain
 - Muscle weakness
 - Muscle soreness

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS308460	OPPO 2321 ELASTIC KNEE SUPPORT - SML
FBS308450	OPPO 2321 ELASTIC KNEE SUPPORT - MED
FBS308440	OPPO 2321 ELASTIC KNEE SUPPORT - LGE
FBS308470	OPPO 2321 ELASTIC KNEE SUPPORT - XL
FBS308480	OPPO 2321 ELASTIC KNEE SUPPORT - XXL