

Aidacare Isolation Cart



- The Aidacare Isolation Cart offers all the desirable features for high-risk infection control areas
- Aidacare carts are designed to offer an attractive, practical and durable cart with multiple configurations to enable the efficient organisation and delivery of medication and equipment
- Smooth, brightly coloured exterior
- Configurable large-volume deep drawers
- Capable of supporting heavier loads with its strong robust construction, this will be your cart of choice
- Easy roll 125 mm castors
- **Standard Inclusions:** waste bin, extended surface and dividers for all drawers
- This value packed cart offers a comprehensive accessory list
- Standard drawer configuration detailed below
- **Small Isolation Cart:**
 - 1 medium 6' drawer (16 cm)
 - 2 large 9' drawers (25 cm)
- **Medium Isolation Cart:**
 - 3 large 9' drawers (25 cm)



Accessories

Click on the cart icon to display the range of optional accessories

Warranty Information

- 5 Years Warranty
- Accessories: 1 Year Warranty

Code	Display Name	Total Length	Total Width	Total Height	Colour
WAT909620	Aidacare Isolation Cart - Yellow - Medium - 100cm / 9H	790 mm	610 mm	1000 mm	Yellow
WAT909610	Aidacare Isolation Cart - Yellow - Small - 90cm / 8H	790 mm	610 mm	900 mm	Yellow