

Liko Golvo 9000 Lifter



- Vertical lifting motion with flexible lift strap is ideal for horizontal lifting, ambulation and lifting from the floor
- Parallel base-width adjustment contributes to easier lifting of patients or residents to/from chairs, beds and the floor
- Armrests help to support patients during ambulation, as well as making it easier for caregivers to manoeuvre the lift
- User display highlights device usage, including lift count, service required indicator and overloads
- Optional quick-release hook system means caregivers can easily and rapidly change a sling bar or install a scale
- Ergonomic hand control for lift operation displays both battery usage and overload indicators
- An extensive assortment of accessories is available including slings, slingbars and scales

Warranty Information

- Frames: 10 Years
- Electrical Parts: 2 Years

 Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material
LSS393100	Liko Golvo 9000 Lifter	1185 mm	735 - 1020 mm	1455 - 2090 mm	200 kg	Aluminium
LSS393150	Liko Golvo 9000 - Leg Protector Set -	-	-	-	-	-