


Liko Golvo 9000 Lifter



- Vertical lifting motion with flexible lift strap is ideal for horizontal lifting, ambulation and lifting from the floor
- Parallel base-width adjustment contributes to easier lifting of patients or residents to/from chairs, beds and the floor
- Armrests help to support patients during ambulation, as well as making it easier for caregivers to manoeuvre the lift
- User display highlights device usage, including lift count, service required indicator and overloads
- Optional quick-release hook system means caregivers can easily and rapidly change a sling bar or install a scale
- Ergonomic hand control for lift operation displays both battery usage and overload indicators
- An extensive assortment of accessories is available including slings, slingbars and scales

 Contract approved

ndis

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Warranty Information

- Frames: 10 Years
- Electrical Parts: 2 Years

| Code | Display Name | Total Length | Total Width | Total Height | Safe Working Load | Material |
|-----------|---------------------------------------|--------------|---------------|----------------|-------------------|-----------|
| LSS393100 | Liko Golvo 9000 Lifter | 1185 mm | 735 - 1020 mm | 1455 - 2090 mm | 200 kg | Aluminium |
| LSS393150 | Liko Golvo 9000 - Leg Protector Set - | - | - | - | - | - |