$\bigcirc$ 



## **TabTimer 8 Alarm Timer**



## **TabTimer 8 Alarm Timer**

Set up to 8 alarms per day. Featuring either an aural or visual alert, magnetised clip, and countdown timer.

- Place on the fridge or other convenient location for medicine or other reminder
- Up-to 8 alarms per day Easy to set and use
- No need to reprogramme, the set alarms will repeat daily
- Countdown timer for regular doses (up to 99 hours and 59 mins)
- Alarm sounds for 60 seconds when reminder time arrives
- 3 large buttons control all functions
- Large easy-read LCD
- 12 hour clock with 'PM' indicator for PM or 24 hour display function (military style)
- Magnet on back, stand on desk, or hang on wall
- 1 x AAA battery included
- Unit Dimensions: 80 x 56 x 20 mm
- Volume approx. 40-45 db
- Weight 59 g
- 12 Month Limited Warranty

Code

**Display Name** 

DLG282000 TabTimer 8 Alarm Timer