

Belinda Bedside Cabinet



- The Belinda Bedside Cabinet delivers simple, understated styling and ample storage
- Three drawers help to reduce clutter and keep essentials beside the bed.
- Deep bottom drawer is ideal for larger items.
- Raised safety edging on the cabinet top prevents items from sliding off.
- Rear towel rail offers a practical and tidy solution for compact living.
- Slimline handles and smooth drawer action offer function and practicality.
- A versatile bedside cabinet to suit most settings and decor.
- Laminate finish is durable and easy to care for.
- Select from our range of laminate colour options.
- Optional hutch provides added shelf space along with practical storage for books and magazines.



Code	Display Name	Total Length	Total Width	Total Height
BEA826310	Belinda Bedside Cabinet - Laminate	450 mm	465 mm	810 mm
BEA826315	Belinda Bedside Cabinet w/Hutch - Laminate	450 mm	465 mm	1020 mm