

## INVISA-BEAM Dual Beam Bed Monitor



- The INVISA-BEAM Dual Beam Bed Monitor is an ideal aid for preventing falls from a bed
- Monitors both sides of the bed with precise invisible beams
- Alerts the nurse at the earliest stage, when the patient attempts to leave the bed
- Assists in patient security by alerting the staff of any unauthorised entry to the bed
- Watch [informational video here](#)
- With a height adjustable stand, it can be used for normal and lowered bed
- Width adjustable to cater for single and king single beds
- As it is a beam that monitors exit or entry to the bed, it can be utilised for bariatric or low weight clients
- Safe and effective for clients positioned on alternating mattresses
- Invisa-beam<sup>®</sup> uses heavy duty, industrial German sensors, which provide long trouble-free life
- Can be utilised with any of the Invisa-Beam wireless alarms as well as connected to nurse call systems
- Optional Wall bracket allows for secure attachment to wall and easy transfer of unit from room to room

 Contract approved

### Support @ Home

**Function:** Signalling, alarming and localisation

**Category:** Prescribed

Click [here](#) for more information

## Warranty Information

- 12 Months
- 24 Months Electronics

<b>Code</b>	<b>Display Name</b>
BEA007216	Dual Beam Monitor w/ Mono & Stereo Plug
BEA007217	Dual Beam Monitor w/ Remote Alarm
BEA007218	Dual Beam Monitor w/ INVISALERT