

McDavid PSII Knee Hinged Thermal Slip-on Support



- Delivers firm compression and soft tissue support
- Made of medical grade lightweight neoprene that retains body warmth
- Contoured design for a better fit with two medial and two lateral medical steel stays
- Top and bottom anchor straps to minimise slippage

✓ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

| Code | Display Name | Knee Circumference |
|------------|--|--------------------|
| FBS3015295 | McDavid - PSII - Knee Hinged - Thermal Slip-on - O/Patella - XXL | 510 - 560 mm |