

## Nosey Cup



- These cups are designed for added comfort and convenience
- A special cutout helps to maintain the proper head and neck positioning when swallowing, making drinking even easier
- Weight: 45 g

### Warranty Information

- 1 Year



Contract approved

**ndis**

**Item Code:**

05\_150900319\_0103\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

| Code      | Display Name                | Capacity |
|-----------|-----------------------------|----------|
| DLK275480 | Nosey Cup 115ml - Blue Tint | 115 ml   |
| DLK275490 | Nosey Cup 225ml - Blue Tint | 225 ml   |