

## Nosey Cup

- These cups are designed for added comfort and convenience
- A special cutout helps to maintain the proper head and neck positioning when swallowing, making drinking even easier
- Weight: 45 g

 Contract approved

**ndis**

**Item Code:**

05\_150900319\_0103\_1\_2

Click [here](#) for more information



### Warranty Information

- 1 Year

### Compliance

No applicable compliance standard

Code	Display Name	Unit Weight	Capacity
DLK275480	Nosey Cup 115ml - Blue Tint	45 g	115 ml
DLK275490	Nosey Cup 225ml - Blue Tint	45 g	225 ml