

## Aspire Deluxe General Purpose Sling with Head Support - Polyester



- Australian standards tested and registered with the TGA
- Versatile sling offering body contoured fabric shaping for maximum patient comfort and support
- This sling has a comfortable padded head support with rigidisers for posture management
- Quality - manufactured in the UK using European materials, computer guided stitching and CNC fabric cutting
- Durable - loop webbing has a breaking strain in excess of one tonne and is designed so that it will not tangle or “cork-screw” during cleaning
- User Comfort & Dignity - features additional padding in weight bearing areas and option to position strap around the outside of legs for more discreet positioning (applicable models)
- Serial Numbered - slings are embroidered with a unique serial number to maintain traceability and compliance at all times
- Service & Maintenance - service record labels ensure you can quickly and easily identify service history
- Easy Sling Identification - all slings are colour coded by size and have reinforced edging for longevity
- Easy Clean - machine washable up to 85°C temperature
- For FIM score 1 or 2

### Warranty Information

- 1 Year

- ✓ Contracted DVA product
- ✓ Contract approved

### ndis

#### Item Code:

05\_123603164\_0105\_1\_2

#### Sling – Composite:

05\_421\_0105\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS ISO 10535:2011

| Code       | Display Name   | Safe Working Load | Material          |
|------------|--|-------------------|-------------------|
| LSS390740  | Deluxe General Purpose Sling with Head Support - Polyester - Small       | 200 kg            | Polyester         |
| LSS390745  | Deluxe General Purpose Sling with Head Support - Polyester - Medium      | 200 kg            | Polyester         |
| LSS390750  | Deluxe General Purpose Sling with Head Support - Polyester - Large       | 300 kg            | Polyester         |
| LSS390755  | Deluxe General Purpose Sling with Head Support - Polyester - Extra Large | 300 kg            | Polyester         |
| LSS3907551 | Deluxe General Purpose Sling with Head Support - Polyester - XXL         | 300 kg            | Polyester (Loops) |
| LSS3907552 | Deluxe General Purpose Sling with Head Support - Polyester - XXXL        | 300 kg            | Polyester (Loops) |