

## Aspire Deluxe General Purpose Sling with Head Support - Polyester



- Australian standards tested and registered with the TGA
- Versatile sling offering body contoured fabric shaping for maximum patient comfort and support
- This sling has a comfortable padded head support with rigidisers for posture management
- Quality - manufactured in the UK using European materials, computer guided stitching and CNC fabric cutting
- Durable - loop webbing has a breaking strain in excess of one tonne and is designed so that it will not tangle or “cork-screw” during cleaning
- User Comfort & Dignity - features additional padding in weight bearing areas and option to position strap around the outside of legs for more discreet positioning (applicable models)
- Serial Numbered - slings are embroidered with a unique serial number to maintain traceability and compliance at all times
- Service & Maintenance - service record labels ensure you can quickly and easily identify service history
- Easy Sling Identification - all slings are colour coded by size and have reinforced edging for longevity
- Easy Clean - machine washable up to 85°C temperature
- For FIM score 1 or 2

### Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Lifting persons

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_123603164\_0105\_1\_2

**Sling – Composite:**

05\_421\_0105\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS390740	Deluxe General Purpose Sling with Head Support - Polyester - Small	200 kg	Polyester
LSS390745	Deluxe General Purpose Sling with Head Support - Polyester - Medium	200 kg	Polyester
LSS390750	Deluxe General Purpose Sling with Head Support - Polyester - Large	300 kg	Polyester
LSS390755	Deluxe General Purpose Sling with Head Support - Polyester - Extra Large	300 kg	Polyester
LSS3907551	Deluxe General Purpose Sling with Head Support - Polyester - XXL	300 kg	Polyester (Loops)
LSS3907552	Deluxe General Purpose Sling with Head Support - Polyester - XXXL	300 kg	Polyester (Loops)