

Jobst Compression Stockings - Men



- Jobst men compression stockings combine style with effective compression support
- Attractive ribbed design with all-day non-constricting comfort grip band keeps socks up without binding or pinching
- Special soft micro fibres wick away moisture six times faster than standard nylon, keeping feet cooler and dryer longer
- Special formula fights odour by inhibiting growth of fungi and bacteria on the socks
- Durable, reinforced heel resists shoe abrasion, strong yarns resist wear and tear

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Warranty Information

- 1 Year

| Code | Display Name | Colour | Calf/Arm Circumference | Knee Circumference | Ankle Circumference |
|-----------|--|--------|------------------------|--------------------|---------------------|
| PTS640040 | Jobst Compression Stockings - Knee - 15-20mmHg - C/Toe - Sml | Black | 28 - 38 cm | 40 -62 cm | 18 -21 cm |
| PTS640530 | Jobst for Men - Knee - 20-30mmHg - C/Toe - Sml | Black | 28-38cm | 40-62cm | 18-21cm |
| PTS640540 | Jobst for Men - Knee - 20-30mmHg - C/Toe - Med | Black | 30-42cm | 46-70cm | 21-25cm |
| PTS640060 | Jobst Compression Stockings - Knee - 15-20mmHg - C/Toe - Lge | Black | 32 - 46 cm | 54 - 78 cm | 25 - 29 cm |
| PTS640550 | Jobst for Men - Knee - 20-30mmHg - C/Toe - Lge | Black | 32-46cm | 54-78cm | 25-29cm |
| PTS640560 | Jobst for Men - Knee - 20-30mmHg - C/Toe - XL | Black | 34-50cm | 60-81cm | 29-33cm |

| Code | Display Name | Colour | Calf/Arm Circumference | Knee Circumference | Ankle Circumference |
|-----------|---|--------|------------------------|--------------------|---------------------|
| PTS640830 | Jobst for Men - Knee - 30-40mmHg - C/Toe - Sml | Black | 28-38cm | 40-62cm | 18-21cm |
| PTS640840 | Jobst for Men - Knee - 30-40mmHg - C/Toe - Med | Black | 30-42cm | 46-70cm | 21-25cm |
| PTS640850 | Jobst for Men - Knee - 30-40mmHg - C/Toe - Lge | Black | 32-46cm | 54-78cm | 25-29cm |
| PTS640860 | Jobst for Men - Knee - 30-40mmHg - C/Toe - XL | Black | 34-50cm | 60-81cm | 29-33cm |
| PTS640601 | Jobst for Men - Thigh - 20-30mmHg - C/Toe - Sml | Black | 28-38cm | 40-62cm | 18-21cm |
| PTS640602 | Jobst for Men - Thigh - 20-30mmHg - C/Toe - Med | Black | 30-42cm | 46-70cm | 21-25cm |
| PTS640604 | Jobst for Men - Thigh - 20-30mmHg - C/Toe - XL | Black | 34-50cm | 60-81cm | 29-33cm |
| PTS640605 | Jobst for Men - Thigh - 30-40mmHg - C/Toe - Sml | Black | 28-38cm | 40-62cm | 18-21cm |
| PTS640606 | Jobst for Men - Thigh - 30-40mmHg - C/Toe - Med | Black | 30-42cm | 46-70cm | 21-25cm |
| PTS640607 | Jobst for Men - Thigh - 30-40mmHg - C/Toe - Lge | Black | 32-46cm | 54-78cm | 25-29cm |
| PTS640608 | Jobst for Men - Thigh - 30-40mmHg - C/Toe - XL | Black | 34-50cm | 60-81cm | 29-33cm |