

Aidacare Aspire A150F Aluminium Folding Lifter



The Aspire 150F is a lightweight and portable folding patient lifter designed for home and community environments.

- **Folds for Transport and Storage:** Folds quickly and conveniently for storage or transport in the boot of the car
- **Ergonomic Handset:** The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- **Off Set Leg Design:** Provides additional accessibility around chairs and furniture in the home
- **Leg Function:** Simple foot operation to open or close legs for lifting stability or accessibility
- **Quick Release Yoke Attachment:** One touch release and attach mechanism, designed for the busy caregiver. Padded yoke for added patient safety
- **Build Quality:** Strong powder coated aluminium frame with industry leading Linak electronics and warranty provision
- **Safety:** Critical features include emergency stop button, emergency lowering and antientrapment technology, with both patient and caregiver safety at heart
- **Intelligent Diagnostics:** State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- **Tente Casters:** 2 x 75 mm front castors for maximum under bed clearance and accessibility. 2 x 100 mm rear locking castors for safety and manoeuvrability
- Item Weight: 32.9 kg
- Tested and compliant to AS/NZS ISO 10535:2011 (the Australian standard for hoists)

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_123604167_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Accessories

Click on the cart icon to display the range of optional accessories

Prescribing Therapists: Click here to download our
[Clinical Assistance Guide](#)

Warranty Information

- Frame: 5 Years
- Electronics: 3 Years
- Battery and Handset: 1 Year

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Safe Working Load	Material
LSS390600	Aidacare Aspire A150F Folding Patient Lifter	4559	1	1220 - 1280 mm	660 - 1040 mm	1873 mm	428 - 1640 mm	580 - 950 mm	150 kg	Aluminium