

Aidacare Aspire A150F Aluminium Folding Lifter



The Aspire 150F is a lightweight and portable folding patient lifter designed for home and community environments.

- **Folds for Transport and Storage**: Folds quickly and conveniently for storage or transport in the boot of the car
- **Ergonomic Handset**: The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- Off Set Leg Design: Provides additional accessibility around chairs and furniture in the home
- **Leg Function**: Simple foot operation to open or close legs for lifting stability or accessibility
- Quick Release Yoke Attachment: One touch release and attach mechanism, designed for the busy caregiver. Padded yoke for added patient safety
- Build Quality: Strong powder coated aluminium frame with industry leading Linak electronics and warranty provision
- Safety: Critical features include emergency stop button, emergency lowering and antientrapment technology, with both patient and caregiver safety at heart
- Intelligent Diagnostics: State of the art Linak Jumbo Care control system.
 Tracks usage and service requirements for safety and compliance
- Tente Casters: 2 x 75 mm front castors for maximum under bed clearance and accessibility. 2 x 100 mm rear locking castors for safety and manoeuvrability
- Item Weight: 32.9 kg
- Tested and compliant to AS/NZS ISO 10535:2011 (the Australian standard for hoists)

Ocontracted DVA product

⊘Contract approved

ndis

Item Code:

05_123604167_0105_1_2

Hoist -

Composite:05_420_0105_1_2 Click **here** for more information

Compliance

AS/NZS ISO 10535:2011

Accessories

Click on the cart icon to display the range of optional accessories



Prescribing Therapists: Click here to download our

Clinical Assistance Guide

Warranty Information

Frame: 5 YearsElectronics: 3 Years

• Battery and Handset: 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Safe Working Load	Material
LSS390600 Aida	acare Aspire A150F Folding Patient Lif	ter 1220 - 1280 mm	660 - 1040 mm	1873 mm	428 - 1640 mm	580 - 950 mm	150 kg	Aluminium