

## Aidacare Aspire A150F Aluminium Folding Lifter



The Aspire 150F is a lightweight and portable folding patient lifter designed for home and community environments.

- **Folds for Transport and Storage:** Folds quickly and conveniently for storage or transport in the boot of the car
- **Ergonomic Handset:** The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- **Off Set Leg Design:** Provides additional accessibility around chairs and furniture in the home
- **Leg Function:** Simple foot operation to open or close legs for lifting stability or accessibility
- **Quick Release Yoke Attachment:** One touch release and attach mechanism, designed for the busy caregiver. Padded yoke for added patient safety
- **Build Quality:** Strong powder coated aluminium frame with industry leading Linak electronics and warranty provision
- **Safety:** Critical features include emergency stop button, emergency lowering and antientrapment technology, with both patient and caregiver safety at heart
- **Intelligent Diagnostics:** State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- **Tente Casters:** 2 x 75 mm front castors for maximum under bed clearance and accessibility. 2 x 100 mm rear locking castors for safety and manoeuvrability
- Item Weight: 32.9 kg
- Tested and compliant to AS/NZS ISO 10535:2011 (the Australian standard for hoists)

✓ Contracted DVA product

✓ Contract approved

### Support @ Home

**Function:** Lifting persons

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_123604167\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS ISO 10535:2011

### Accessories

Click on the cart icon to display the range of optional accessories

**Prescribing Therapists:** Click here to download our  
[Clinical Assistance Guide](#)

## Warranty Information

- Frame: 5 Years
- Electronics: 3 Years
- Battery and Handset: 1 Year

| Code      | Display Name                                 | Total Length   | Total Width   | Total Height | Lift Range    | Leg Gait Range | Safe Working Load | Material  |
|-----------|--|----------------|---------------|--------------|---------------|----------------|-------------------|-----------|
| LSS390600 | Aidacare Aspire A150F Folding Patient Lifter | 1220 - 1280 mm | 660 - 1040 mm | 1873 mm      | 428 - 1640 mm | 580 - 950 mm   | 150 kg            | Aluminium |