

iCare Over Bed Pole



- Pull up bar to help with turning, patient transfer, and bed care. A stable design with heavy-duty steel, compatible with all ICARE Beds.
- Can be positioned anywhere on side or end of bed
- Free standing design
- Strong bracing brackets
- Adjustable handle height
- Anti-slip strap cradle

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material
BEA016420	iCare Over Bed Pole	900 mm	750 mm	1900 mm	100 kg	Steel