

## DONJOY DP151EB TRIZONE ELBOW SLEEVE



- Three compression zones deliver peak performance.
- Silicone banding injected directly over the sleeve.
- Targets pain from golf and tennis elbow.
- Stabilizes the joint and supports the bicep and triceps muscles.
- All-natural Carbonized Bamboo is thermal regulating and anti-bacterial to help eliminate odours.
- Anti-migration technology helps to keep the elbow sleeve secured to the bicep and eliminates slippage.
- Ergonomic fit.

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Upper limb orthoses

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name
FBS306100	DONJOY DP151EB TRIZONE ELBOW SLEEVE - SML
FBS306090	DONJOY DP151EB TRIZONE ELBOW SLEEVE - MED
FBS306080	DONJOY DP151EB TRIZONE ELBOW SLEEVE - LGE
FBS306110	DONJOY DP151EB TRIZONE ELBOW SLEEVE - XL