

Over Toilet Aid - Maxi



- Extra wide and reinforced mild steel tubing frame accomodating maximum user weights up to 300 kg
- Fixed height frame for further strength and stability
- Removable lift back seat for easy cleaning, comes standard with splash guard
- Fitted with non-slip rubber tips
- Zinc treated anti corrosive frame

 Contract approved

Support @ Home

Function: Toileting

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 1 Year

Compliance

AS/NZS 3973:2009

Code	Display Name	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height
BTT146500	Over Toilet Aid - Maxi 600mm Wide	300 kg	5.5 kg	420 mm	600 mm	520 mm
BTT146600	Over Toilet Aid - Maxi 650mm Wide	300 kg	6.5 kg	420 mm	650 mm	520 mm