

## Over Toilet Aid - Maxi



- Extra wide and reinforced mild steel tubing frame accomodating maximum user weights up to 300 kg
- Fixed height frame for further strength and stability
- Removable lift back seat for easy cleaning, comes standard with splash guard
- Fitted with non-slip rubber tips
- Zinc treated anti corrosive frame

### Warranty Information

- 1 Year



Contract approved

### Support @ Home

**Function:** Toileting

**Category:** Prescribed

Click [here](#) for more information

### Compliance

AS/NZS 3973:2009

Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
BTT146500	Over Toilet Aid - Maxi 600mm Wide	420 mm	600 mm	520 mm	300 kg	5.5 kg
BTT146600	Over Toilet Aid - Maxi 650mm Wide	420 mm	650 mm	520 mm	300 kg	6.5 kg