

Limbkeepers - Leg Sleeve



- Non-compression soft & breathable knit
- Latex free
- UPF 50+ Sun Protection
- Anti-microbial - silver infused yarns help curb odour causing bacteria (sleeves only)
- Wicking - channels moisture away from skin
- Seamless - no fasteners
- Stays in place - no adhesives
- Form fitting with 4-way rebound stretch
- Unique blend of performance yarns - keep cool when it's warm/warm when it's cool
- Non-restricting
- Interchangeable design for a perfect fit
- Colourfast - easy care/quick dry
- Soft and breathable knit
- One size fits most
- Unisex
- Sold as Single Indications:
- Designed to be used as under sleeve for knee brace, protection for post knee surgery scarring and healing, or from hazards of daily activities

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Body-worn assistive products for body protection

Category: Under advice

Click [here](#) for more information

Code	Display Name	Colour
PTA537560	Limbkeepers - Limb Protection - Leg Sleeve - Heavy(Regular)	Standard
PTA537570	Limbkeepers - Limb Protection - Leg Sleeve - Heavy(Regular) - Black	Black

Code	Display Name	Colour
PTA537580	Limbkeepers - Limb Protection - Leg Sleeve - Heavy(Regular) - White	White
PTA537610	Limbkeepers - Limb Protection - Leg Sleeve - Light	Standard
PTA537620	Limbkeepers - Limb Protection - Leg Sleeve - Light - Black	Black
PTA537630	Limbkeepers - Limb Protection - Leg Sleeve - Light - White	White