

Oxford Presence Patient Lifter



- Mobile patient lifter for transfer to a seated or lying position
- High lifting range accommodates transfers to higher surfaces
- Large load capacity of 227 kg
- Robust epoxy powder coated frame construction makes it ideal for facilities
- Electric hand control operation
- Emergency lowering and raising
- Over-sized push handles for excellent manoeuvrability
- 100 mm rear locking castors and 100 mm front twin swivel castors
- Electric leg gating
- Suits a large range of patient slings
- Supplied with battery and charger
- Pictured with optional 4 point pivot frame (LSS377395)
- Item Weight: 40 kg



Contract approved

ndis

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 2 Years

| Code | Display Name | Total Length | Total Width | Total Height | Lift Range | Leg Gait Range | Electric Leg | Safe Working Load | Material |
|-----------|-------------------------|--------------|-------------|----------------|---------------|----------------|--------------|-------------------|-----------|
| LSS377325 | Presence Patient Lifter | 1460 mm | 670 mm | 1460 - 2061 mm | 440 - 1910 mm | 550 - 1020 mm | True | 227 kg | Aluminium |