

ROHO Standard Cushion Cover



permobil

ROHO

- The ROHO® Standard Series Cushion Cover is a textile enclosure, intended to protect the cushion
- Two -way stretch nylon/spandex top, polyester spacer fabric sides with zipper, and a polyester/rayon non-skid PVC coated bottom with hook and loop fasteners to help secure the cushion in the chair

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Total Width
PTC596200	ROHO Cushion Cover - HP - Standard - 18.5x18.5in	-	-
PTC596300	ROHO Cushion Cover - LP - Standard - 18.5x18.5in	-	-
PTC596400	ROHO Cushion Cover - HP - Standard - 18.5x20in	-	-
PTC596500	ROHO Cushion Cover - HP - Standard - 18.5x22in	-	-
PTC596600	ROHO Cushion Cover - HP - Standard - 20x18.5in	-	-
PTC596700	ROHO Cushion Cover - HP - Standard - 20x20in	-	-
PTC596720	ROHO Cushion Cover - LP - Standard - 20.5x18.5in	-	-

Code	Display Name	Total Length	Total Width
PTC596750	ROHO Cushion Cover - HP - Standard - 22x18.5in	-	-
PTC596800	ROHO Cushion Cover - HP - Standard - 22x20in	-	-
PTC596900	ROHO Cushion Cover - HP - Standard - 11x13in	-	-
PTC597000	ROHO Cushion Cover - HP - Standard - 13.25x16.75in	-	-
PTC597100	ROHO Cushion Cover - HP - Standard - 15x15in	-	-
PTC597150	ROHO Cushion Cover - LP - Standard - 15x15in	-	-
PTC597200	ROHO Cushion Cover - HP - Standard - 15x16.75in	-	-
PTC597300	ROHO Cushion Cover - LP - Standard - 15x16.75in	-	-
PTC597350	ROHO Cushion Cover - HP - Standard - 15x18.5in	-	-
PTC597400	ROHO Cushion Cover - HP - Standard - 16.75x18.5in	-	-
PTC597500	ROHO Cushion Cover - LP - Standard - 16.75x18.5in	-	-
PTC597550	ROHO Cushion Cover - HP - Standard - 16.75x20in	-	-
PTC597600	ROHO Cushion Cover - LP - Standard - 16.75x20in	-	-
PTC597850	ROHO Cushion Cover - HP - Standard - 16.75x15in	-	-
PTC597900	ROHO Cushion Cover - LP - Standard - 16.75x16.75in	-	-
PTA523700	Custom Made ROHO Cover (As per script)	-	-
PTC555740	ROHO Cover Waterproof - LP - 20inx18.5in	18.5 in	20 in

