


Pride 660 Powerlift Recliner - Dual Motor



- To suit the bariatric user
- Engineered furniture laminate/hardwood frame
- Patented, quiet & smooth lift system
- Large side pocket
- Hand control can be easily moved to left-hand side
- Integrated emergency battery backup (batteries not included)
- External transformer for easy access to battery backup

Warranty Information

- 2 Years

 Contract approved

ndis

Item Code:

05_180903180_0103_1_2

Click [here](#) for more information

Compliance

ISO13934-1 / ISO13937-1 / BS6807
/ BS7175

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

| Code | Display Name | RRP | Box Quantity | Seat Depth | Seat Width | Backrest Height | Seat To Floor Height | User Height | Safe Working Load | Colour |
|-----------|----------------------|------|--------------|------------|------------|-----------------|----------------------|-------------|-------------------|--------------|
| CHP193702 | 660 Recliner - Latte | 3029 | 1 | 510 mm | 510 mm | 635 mm | 480 mm | < 1630 mm | 170 kg | Latte Fabric |