

Liko LiftPants (Mod. 92/920)



- Enables gait training and rehab to occur, whilst offering patients an increased sense of security
- Offer freedom of movement while relieving the burden of body weight
- Provide necessary support to allow patient to move on his or her own with confidence
- Can be used on overhead or mobile patient lifters
- Available in Mesh



ndis

Item Code:

05_123603166_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Warranty Information

- 1 Year

Compliance

AS/NZS ISO 10535:2011



Code	Display Name	Safe Working Load	Material
LSS394340	Lift Pants Mod 92 Walking Sling - Mesh - SMALL	200 kg	Mesh (Loops)
LSS394350	Lift Pants Mod 92 Walking Sling - Mesh - MEDIUM	200 kg	Mesh (Loops)
LSS394360	Lift Pants Mod 92 Walking Sling - Mesh - LARGE	200 kg	Mesh (Loops)
LSS394370	Lift Pants Mod 92 Walking Sling - Mesh - XLARGE	500 kg	Mesh (Loops)
LSS394380	Lift Pants Mod 92 Walking Sling - Mesh - XXLARGE	500 kg	Mesh (Loops)
LSS394390	Lift Pants Mod 92 Walking Sling - Mesh - XXXLARGE	500 kg	Mesh (Loops)