

Liko Amputee Sling (Mod. 75)



- Designed for safe and secure lifting of patients with amputated legs
- Supports patient with single and double leg amputations as well as patients with contractures
- Can also be used for non-amputee patients who are slender and require a smaller sling
- Offers a sitting posture that is slightly reclined with leg supports applied parallel to each under both thighs
- Model 75 offers head support and is available in Polyester

Warranty Information

- 1 Year



 Contract approved

ndis

Item Code:

05_123603166_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

| Code | Display Name | Safe Working Load | Material |
|-----------|--|-------------------|-------------------|
| LSS394520 | Amputee Highback Sling - Reinforced Leg Support - Polyester - SMALL | 200 kg | Polyester (Loops) |
| LSS394530 | Amputee Highback Sling - Reinforced Leg Support - Polyester - MEDIUM | 200 kg | Polyester (Loops) |
| LSS394540 | Amputee Highback Sling - Reinforced leg Support - Polyester - LRG | 200 kg | Polyester (Loops) |
| LSS394550 | Amputee Highback Sling - Reinforced Leg Support - Polyester - XLARGE | 300 kg | Polyester (Loops) |