

Lopital Shower Commode - Reflex



- Electrically operated height-adjustable settings (variable)
- Electrically operated tilt mechanism
- Easy-to-use touch controls on backrest
- Height-adjustable headrest
- Ergonomically designed backrest and seat
- Height-adjustable legrest and calf protection
- Foldable swivel armrests with front locking
- Suitable for use in the shower and toilet
- Can be moved into position over toilet
- Slide for stainless steel commode pan
- Four plastic castors with brakes
- Easy to recharge with battery charger provided



Warranty Information

- 2 Years

| Code | Display Name | Safe Working Load | Unit Weight | Seat Depth | Seat Width | Seat To Floor Height | Seat To Footrest |
|-----------|-------------------------|-------------------|-------------|------------|------------|----------------------|------------------|
| BTC092400 | Shower Commode - Reflex | 135 kg | 55 kg | 48.5 cm | 48 cm | 49 - 101 cm | 49 - 101 cm |