

Isotoner - Therapeutic pressure gloves



- Washable lycra and spandex blend gives Isotoner gloves a comfortable fit while controlling swelling
- Use this glove anytime day or night and the gentle, even compression will soothe and reduce swollen aching hands
- Seamless support - seams are on the outside and away from pain-sensitive areas like the base of the thumb
- 2-way stretch provides gentle compression
- Isotoner therapeutic glove compression is considered moderate: 23 - 32mm Hg

 Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name
FBS301159	Isotoner - Therapeutic pressure Open Finger - Sml
FBS301160	Isotoner - Therapeutic pressure Open Finger - Med
FBS301161	Isotoner - Therapeutic pressure Open Finger - Lge