

Male Insert Pads



Conni.

- Soft micro-fibre lining rapidly absorbs and feels dryer, flexible PUL waterproof is cool and comfortable with an elastic clip strap for secure application
- Specifically shaped for men, supplied in a twin pack with a convenient waterproof zip bag
- Designed for light incontinence (up to 250 mls over 4 hrs)

Warranty Information

- 1 Year

 Contract approved

Support @ Home

Function: Absorbing products to contain urine and faeces

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

03_093021078_0103_1_1

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Total Length
CON228570	Insert Pad Mens Washable Twin Pack w/Bag	390 mm