

Aspire Bath Transfer Bench - Plastic Seat



- The Aspire Transfer Bench has a textured moulded plastic seat & backrest
- Fully height adjustable with a side arm for secure transfers
- Adjustable backrest to suit right or left facing tubs
- Suction cup foot pieces secure the bench to the tub floor
- Lightweight aluminium frame construction
- **Aspire Extension Legs - [Set of 4](#)**

Each set includes:

- 2 x suction tips (for use on smooth surfaces)
- 2 x standard tips (for general use)

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 1 Year



Code	Display Name	Safe Working Load	Material	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
BTS108200	Aspire Bath Transfer Bench - Plastic Seat	181 kg	Plastic Seat Aluminium Frame	4.55kg	405 mm	670 mm	455-570 mm	455-570 mm
BTS108332	Extension Legs - 100mm - w/2 Suction & 2 Std Tips - Suits BTS108200	-	-	-	-	-	-	-