

Rolyan AirPro Air Graduate Splint



- Air Therapy-low load prolonged stretch
- Conforms to the shape of the hand accommodating different degrees of finger contracture
- Ideal for severely contracted fingers
- Removable cone shaped palmar air bladder can be used independently of wrist support
- Extra cover included
- Universal size

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name
FBS301494L	Rolyan AirPro Air Graduate Splint - Left
FBS301494	Rolyan AirPro Air Graduate Splint - Right