

Aspire Mobility Scooter Rear Bag



- The Aspire Scooter rear bag is a convenient storage location to hold shopping and personal belongings
- Please check before purchasing if this can be attached to your Aspire Mobility Scooter

✓ Contracted DVA product

✓ Contract approved

Compliance

AS3695.2:2013

Code	Display Name	RRP	Box Quantity	Total Length	Total Height	Seat Depth
SCA668630	Aspire Mobility Scooter Rear Bag	317	1	400 mm	450 mm	250 mm