

BODY ASSIST 388B THERMAL ANKLE WRAP



- General purpose thermal ankle compression support.
- Medical grade CR neoprene is LATEX FREE
- Medical plush inner liner draws perspiration, delivers warmth and provides soft plush comfort against the skin.
- Adjustable closure with contoured design for a snug fit.
- Black looks great as part of any sports uniform. Application:
- Suitable for simple sports strains/sprains.

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Colour
FBS307680	BODY ASSIST 388B THERMAL ANKLE WRAP - BLACK - UNIVERSAL	Black
FBS307690	BODY ASSIST 388S THERMAL ANKLE WRAP - BEIGE - UNIVERSAL	Beige