

## Venapure - Anti Embolic - Compression Stockings



- The VenaPure stockings offer a greater degree of two way combined with an increased end stretch provide for a better anatomical fit to the patient.
- The close knit construction provides higher durability of the material. The material is also latex free
- The inspection toe hole position sits on top of the toe area making foot inspection easier. The change in position of the inspection hole also prevents the toes from protruding from the inspection hole on walking
- Venapure knee stockings are designed to reduce venous stasis below the knee when thigh length styles are medically contraindicated
- Thigh length with grip top - Pressure break at knee, protecting the critical popliteal area and promoting free leg movement in the recumbent patient
- The Venapure thigh high with waist belt is very suitable for thin or obese patients, where stocking slippage is often a problem - One-piece, snap-on belt is easily adjustable
- Available in various sizes, with choice of either short, regular or long fit

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Supporting and stimulating blood circulation

**Category:** Low risk

Click [here](#) for more information

Code	Display Name	Colour
FBS305120	VENAPURE - Below-Knee - Anti-embolic - Toe Hole - Long / Sml - White	White
FBS305130	VENAPURE - Below-Knee - Anti-embolic - Toe Hole - Long / Med - White	White
FBS305140	VENAPURE - Below-Knee - Anti-embolic - Toe Hole - Long / Lge - White	White
FBS305150	VENAPURE - Below-Knee - Anti-embolic - Toe Hole - Long / XL - White	White
FBS305160	VENAPURE - Below-Knee - Anti-embolic - Toe Hole - Long / XXL - White	White

Code	Display Name	Colour
FBS305180	VENAPURE - Below-Knee - Anti-embolic - Toe Hole - Reg / Sml - White	White
FBS305190	VENAPURE - Below-Knee - Anti-embolic - Toe Hole - Reg / Med - White	White
FBS305200	VENAPURE - Below-Knee - Anti-embolic - Toe Hole - Reg / Lge - White	White
FBS305210	VENAPURE - Below-Knee - Anti-embolic - Toe Hole - Reg / XL - White	White
FBS305220	VENAPURE - Below-Knee - Anti-embolic - Toe Hole - Reg / XXL - White	White
FBS305230	VENAPURE - Below-Knee - Anti-embolic - Toe Hole - Reg / XXXL - White	White
PTS643100	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Long / Sml - White	White
PTS643110	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Long / Med - White	White
PTS643120	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Long / Lge - White	White
PTS643130	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Long / XL - White	White
PTS643140	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Long / XXL - White	White
PTS643150	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Reg / Sml - White	White
PTS643160	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Reg / Med - White	White
PTS643170	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Reg / Lge - White	White
PTS643180	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Reg / XL - White	White
PTS643190	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Reg / XXL - White	White
PTS643200	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Short / Sml - White	White

Code	Display Name	Colour
PTS643210	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Short / Med - White	White
PTS643220	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Short / Lge - White	White
PTS643230	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Short / XL - White	White
PTS643240	VENAPURE - Thigh w/ Grip top - Anti-embolic - Toe Hole - Short / XXL - White	White
PTS643250	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - Long / Sml - White	White
PTS643260	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - Long / Med - White	White
PTS643270	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - Long / Lge - White	White
PTS643280	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - Long / XL - White	White
PTS643290	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - Long / XXL - White	White
PTS643300	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - Reg / Sml - White	White
PTS643310	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - Reg / Med - White	White
PTS643320	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - Reg / Lge - White	White
PTS643330	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - Reg / XXL - White	White
PTS643340	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - Reg / XL - White	White
PTS643350	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - SHORT / Sml - White	White
PTS643360	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - SHORT / Med - White	White
PTS643370	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - SHORT / Lge - White	White
PTS643380	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - SHORT /XL - White	White

Code	Display Name	Colour
PTS643390	VENAPURE - Thigh w/ Belt - Anti-embolic - Toe Hole - SHORT / XXL - White	White