

iCare Bed Sheet Set



- Bed sheet designed to fit and move with iCare profiling beds.
- Luxurious 1000 Thread Count
- 100% Cotton Plain Sateen

 Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Code	Display Name
BEA830500	iCare Bed Sheet Set - Queen
BEA830505	iCare Bed Sheet Set - Long Double
BEA830510	iCare Bed Sheet Set - Long Single
BEA830515	iCare Bed Sheet Set - Split King
BEA830520	iCare Bed Sheet Set - Split Queen