

## Maxxcare Heel Pro Evolution



The Maxxcare Heel Pro Evolution provides excellent comfort and pressure reduction for patients with a 'high risk' of developing pressure ulcers to the vulnerable heel area. For patients with pressure damage, the Maxxcare Heel Pro Evolution can be an effective solution for wound healing, as its unique design allows for heel offloading.

- Lightweight 'SmartCell' air technology makes this a practical and easy to wear heel protector.
- Provides excellent levels of comfort and pressure reduction.
- Improved non-slip sole means that you can safely walk short distances while wearing the Maxxcare Heel Pro Evolution.
- Foam free, low friction nylon outer material ensure easy movement over sheets, improving bed mobility and circulation.
- 60 oC (machine washable) without removing the internal air cushions.

✔ Contract approved

### Support @ Home

**Function:** Body-worn assistive products for body protection

**Category:** Under advice

Click [here](#) for more information

Code	Display Name	Ankle Circumference
PTA509350	Maxxcare Heel Pro Evolution - Small	Less than 380 mm
PTA509355	Maxxcare Heel Pro Evolution - Medium	380 - 430 mm
PTA509360	Maxxcare Heel Pro Evolution - XL	430 mm