

SitnStand Lift Assist - Compact



SitnStand Lift Assist Compact is designed for those who have difficulty standing up. This “lift chair in a bag” is Ultra-lightweight and portable and can go wherever you go – parks, restaurants, theaters, arenas, and even on airplanes or mobility scooters.

- Fits any chair, sofa, seat, bench, or scooter outside the home
- [Watch instructional video here.](#)
- Simple to Use: Two-button remote for full control
- VLT™ - Vertical Lift Technology
- Comes with a carry bag for easy transport
- Rechargeable battery goes 4-7 days between recharging
- Ultra-lightweight: only 3kg
- Rise up to 20cm with Ease

Warranty Information

- 2 Year Warranty

✔ Contract approved

Support @ Home

Function: Changing body position

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Safe Working Load	Unit Weight
LSS403690	SitnStand Lift Assist - Compact	121kg	3kg