

## Aidacare AC4 Bed



The Aidacare AC4 bed is the perfect blend of function and durability, providing a reliable sleeping solution for your residents. The gentle gunmetal bed frame coupled with teak bed ends or upholstered surround options, makes it easy to blend with interior styling, delivering a homely ambience. Designed and manufactured to the highest standards with trusted TiMotion electronics for long-term reliability and minimal maintenance.

- Multiple functions and accessories combine to exceed the comfort, support, and safety expectations for a diverse range of users, whilst facility carers benefit from the considered design.
- **Designed for carers and residents** - TiMotion electronics are up to 20% faster than other brands, saving carers time. Simple and intuitive handset is easy to operate.
- **Maintains infection control** - Wipe clean mattress platform allows easy cleaning for infection control compliance.
- **Personalised sleeping solution**
  - Quick attaching and easily removable head and foot boards
  - 200 mm mattress platform extension caters for taller residents\*
  - Optional upholstered Allure surround for an alternative homely style
- **Considered cable management** - Ensures cables are kept safely out of the way, minimising damage and downtime.
- **Ensures wall protection** - Wall locator bar included and bed rises vertically thereby minimising costly wall collisions while maximising functional bedroom space at the foot-end.
- **Enhanced manual handling**
  - 125 mm twin castors: Offer lower rolling resistance for easier bed repositioning
  - Central Locking - Total & Directional: Activated centrally via an easily

✔ Contract approved

### Support @ Home

**Function:** Beds and bed equipment

**Category:** Prescribed

Click [here](#) for more information

### Compliance

AS/NZS IEC 60601-1:2015

IEC 60601-2-52 (Edition 1.1 2015-03)

accessible full-width foot pedal, push down for total lock, or lift up for directional lock to facilitate safe and easy tracking through corridors and doorways.

- Optional Accessories available - Comes standard with Sublime teak laminate head and foot boards
- **Key Points Overview**
  - Premium quality Australian designed beds for Aged Care
  - Designed for carers and residents
  - Enhanced manual handling
  - Multi-functional positioning
  - Durable, reliable, low-maintenance
  - Contemporary design
  - Bespoke finishes available

## Warranty Information

- Frame: 10 Years
- Electronics: 3 Years
- Handset: 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
BEB050400	Aidacare AC4 Bed - Single	2100 mm*	920 mm	300 - 730 mm	220 kg	112 kg
BEB050410	Aidacare AC4 Bed - King Single	2100 mm*	1070 mm	300 - 730 mm	220 kg	116kg