

Rolyan Kwik-Form Functional Resting Splint



- Can be custom fitted without tools or heat - simply bend the wire core to the desired position
- Supports the wrist, thumb and fingers in a functional or resting position
- Features hook and loop closures and removable machine washable cover



Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name
FBS301483	Rolyan Kwik-Form Functional Resting Splint - Right - Sm/Med