

Raised Toilet Seat with Swing Back Arms



- This raised toilet seat features swing back arms and is height adjustable
- Fits most standard toilets
- Fixes onto the toilet bowl for safe transferring (replaces existing seat & lid).
- Plastic moulded armrests
- Height adjustable range: 5, 10 and 15 cm
- Width Between Arms: 410 mm

Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Toileting

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_091218063_0103_1_2

Toilet Aids – Composite:

05_414_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name
------	--------------

BTT150450	Raised Toilet Seat with Swing Back Arms
-----------	---