

Raised Toilet Seat with Swing Back Arms



- **Anti-slip Support Armrests:** Textured armrests allow users to safely grip the armrests when sitting and standing from the Raised Toilet Seat for support.
- **Swing-back for Easy Transfers with Walking Aids:** Swing-back mechanism allows armrest to be moved out of the way to assist transfers from wheelchairs, and walking aids.
- **Hygienic and Easy to Clean:** Including a lid and components made from moulded plastic, it is easy to clean to minimize infection risks.
- **Universal Fit:** Featuring depth and height adjustable mounting bolts, the Aspire Toilet Seat Raiser is versatile to suit almost any toilet.
- **Height Adjustable:** Suitable for users of various heights or with mobility limitations, helping them regain confidence and remain independent at home.

Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Toileting

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_091218063_0103_1_2

Toilet Aids – Composite:

05_414_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name	Total Width	Total Height	Safe Working Load
BTT150450	Raised Toilet Seat with Swing Back Arms	-	-	-
BTT164050	Aspire Toilet Seat Raiser with Swing Back Arms	565mm	370mm	136kg

