

## Cushion - Seat - Memo Gel



- Control your body temperature while seated with this memory gel seat cushion.
- Dual-sided design ensures adaptability for both hot and cool days.
- Allows you to enjoy the cooling effects of the gel-infused memory foam on hot days by placing the silica gel layer on top, or flip it for the comforting, moulding embrace of memory foam on cooler days.
- Incorporates a silica gel layer that actively controls your body temperature.
- Its adaptability means you can bring the cooling, warming and supportive benefits of the cushion with you, making it your perfect travel companion

 Contract approved

### Support @ Home

**Function:** Accessories for wheelchair seating

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Total Height
BEP056214	Memory Gel Seat Cushion	400 mm	400 mm	50 mm