

## Cushion - Seat - Memo Gel



- Control your body temperature while seated with this memory gel seat cushion.
- Dual-sided design ensures adaptability for both hot and cool days.
- Allows you to enjoy the cooling effects of the gel-infused memory foam on hot days by placing the silica gel layer on top, or flip it for the comforting, moulding embrace of memory foam on cooler days.
- Incorporates a silica gel layer that actively controls your body temperature.
- Its adaptability means you can bring the cooling, warming and supportive benefits of the cushion with you, making it your perfect travel companion



Code	Display Name	Total Length	Total Width	Total Height
BEP056214	Memory Gel Seat Cushion	400 mm	400 mm	50 mm