

## BODY ASSIST 409 THERMAL KNEE WRAP



- General purpose thermal knee compression wrap.
- Closed patella.
- Medical grade CR neoprene is LATEX FREE.
- Adjustable closure with contoured design for a snug fit.
- Black looks great as part of any sports uniform.
- Suitable for simple sports strains/sprains.
- Machine washable for easy care.
- Fits left or right knee.
- Retains body warmth & can increase blood circulation. Application:
- To protect the knee during recovery.
- To cushion against minor scrapes and bumps.
- Low level support & protection to painful, tired knees.

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Lower limb orthoses

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name
FBS307710	BODY ASSIST 409 THERMAL KNEE WRAP - UNIVERSAL